Application Handbook Thermal Analysis

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This application handbook presents selected application examples. The experiments were conducted with the utmost care using the instruments specified in the description of each application at METTLER TOLEDO Thermal Analysis Lab in Switzerland. The results have been evaluated according to the current state of our knowledge.

This does not however absolve you from personally testing the suitability of the examples for your own methods, instruments and purposes. Since the transfer and use of an application is beyond our control, we cannot of course accept any responsibility.

When chemicals, solvents and gases are used, general safety rules and the instructions given by the manufacturer or supplier must be observed.

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Thermal analysis is the name given to a group of techniques used to measure the physical and chemical properties of materials as a function of temperature. In all these methods, the sample is subjected to a heating, cooling or isothermal temperature program.

The measurements can be performed in different atmospheres. Usually either an inert atmosphere (nitrogen, argon, helium) or an oxidative atmosphere (air, oxygen) is used. In some cases, the gases are switched from one atmosphere to another during the measurement. Another parameter sometimes selectively varied is the gas pressure.

DSC can also be used in combination with instruments that allow the sample to be simultaneously observed (DSC microscopy) or exposed to light of different wavelengths (photocalorimetry).

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In DSC, the heat flow to and from the sample is measured. DSC can be used to investigate thermal events such as physical transitions (the glass transition, crystallization, melting, and the vaporization of volatile compounds) and chemical reactions. The information obtained characterizes the sample with regard to its thermal behavior and composition. In addition, properties such as the heat capacity, glass transition temperature, melting temperature, heat and extent of reaction can also be determined.

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Conventional DSC employs a linear temperature program. The sample and reference material (or just an empty crucible) are heated or cooled at a linear rate, or in some cases, held at a constant temperature (i. e. isothermally). Often several partial programs or so-called segments are joined together to form a complete temperature program. A typical DSC curve of a polymer is shown schematically in Figure 1.1.

The change in the curve at the beginning of the measurement is due to the initial “startup deflection” (1). In this transient region, the conditions suddenly change from an isothermal mode to a linear heating mode. The magnitude of the startup deflection depends on the heat capacity of the sample and the heating rate.

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At a glass transition (2), the heat capacity of the sample increases and therefore an endothermic step is observed. Cold crystallization(3) occurs above the glass transition. Polymers that crystallize readily are heated to above the melting point and quench-cooled to suppress the formation of crystallites. Such polymers recrystallize above the glass transition. On further heating, melting (4) takes place. At higher temperatures, decomposition (6) begins. Some substances may vaporize(5) between the melting and decomposition.

The type of purge gas used in the experiment often has an influence on the reactions involved.


图 1.1 图示聚合物 DSC 曲线:1 初始启动偏移;2 玻璃化转变;3 冷结晶;4 熔融;5 汽化;6 分解
Fig. 1.1 Schematic DSC curve of a polymer: 1 initial startup deflection; 2 glass transition; 3 cold crystallization; 4 melting; 5 vaporization; 6 decomposition

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Transitions and reactions can be differentiated by cooling the sample and measuring it again-chemical reactions are irreversible whereas crystalline materials melt then crystallize again on cooling or on heating a second time. Glass transitions are also reversible but not the enthalpy relaxation often observed in the first heating measurement of a glass transition.

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Alternating DSC (ADSC) is a particular type of temperature-modulated DSC (TMDSC). In contrast to conventional DSC, the linear temperature program is overlaid with a small periodic temperature change. The temperature program is characterized

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by the underlying heating rate, the temperature amplitude and the duration of the periodically changing temperature (Fig. 1.2). With quasi-isothermal measurements, the underlying heating rate can also be zero.

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Fig. 1.2 Typical ADSC temperature program: β_u is the underlying heating rate, A_T the temperature amplitude, t_p period, The angular frequency ω is defined as $2\pi/P$ where P denotes the period of the sine wave.

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As a result of temperature modulation, the measured heat flow changes periodically. This can be separated into two parts as shown in Figure 1.3. Signal averaging yields the underlying signal (total heat flow), which corresponds to the conventional DSC

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Fig. 1.3 Separation of the measured ADSC curve into the underlying and the periodic signal components

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curve. As additional information, one also obtains the periodic signal component. The reversing heat flow corresponds to the heat flow component that is able to follow the heating rate change directly and is computed from the in-phase heat capacity. The difference between the total heat flow and the reversing heat flow yields the non-reversing heat flow. One advantage of this technique is that it allows processes that occur simultaneously to be separated. For example, the change in heat capacity during a chemical reaction can be measured directly.

The evaluation of the ADSC curves is based on Fourier analysis. The modulus of the complex heat capacity c_p^* is calculated using the equation,

where A_Φ and A_β denote the amplitudes of the modulated heat flow and heating rate, and m the sample mass. The phase angle between the ADSC heat flow signal and the heating rate is used to calculate the in-phase c_p .

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Fig. 1.4 IsoStep temperature program consisting of different isothermal and dynamic segments

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The isothermal segments allow the isothermal drift of the dynamic segments to be corrected. This results in better heat capacity accuracy. The isothermal step may also contain kinetic

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information, for example of a chemical reaction. Heat capacity determinations can be made using a sapphire reference sample, and kinetic effects can be separated from changes in heat capacity.

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TOPEM™ is an advanced temperature-modulated DSC technique that is based on the full mathematical analysis of the response of a DSC (both the apparatus and the sample) to a stochastically modulated underlying temperature program (Fig. 1.5). Due to the randomly distributed temperature pulses, the system is subjected to temperature oscillations over a wide frequency range and not just at one single frequency (ADSC). An analysis of the correlation of the oscillating input signal (heating rate) and the response signal (heat flow) provides much more information than can be obtained using conventional temperature-modulated DSC. Not only can reversing and non-reversing effects be separated, but the quasi-static heat capacity of the sample is also measured and frequency-dependent heat capacity values are determined. This can be used to distinguish between frequency-dependent relaxation effects (e. g. glass transitions) and frequency-independent effects (e. g. chemical reactions) in one single measurement.

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Fig. 1.5. Temperature curve of the furnace set value (black line) in a TOPEM™ method in which the furnace temperature (red curve) generates a heating rate that fluctuates around a mean value.

The heat flow to the sample also fluctuates irregularly as shown in the diagram on the right.

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When a sample is heated, it often begins to lose mass. This loss of mass can result from vaporization or from a chemical reaction in which gaseous products are formed and evolved from the

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sample. If the purge gas atmosphere is not inert, the sample can also react with the gas. In some cases, the sample mass may also increase, e. g. in an oxidation reaction if the product formed is a solid.

In thermogravimetric analysis (TGA), the change in mass of a sample is measured as a function of temperature or time.

TGA provides information on the properties of the sample and its composition. If the sample decomposes as a result of a chemical reaction, the mass of the sample often changes in a stepwise fashion. The temperature at which the step occurs characterizes the stability of the sample material in the atmosphere used.

Figure 1.6 shows a typical TGA curve. The composition of a material can be determined by analyzing the temperatures and the heights of the individual mass steps.

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Fig. 1.6 Schematic TGA curve; 1 loss of mass due to the vaporization of volatile components; 2 pyrolysis in an inert atmosphere; 3 combustion of carbon on switching from an inert to an oxidative atmosphere; 4 residue.

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Volatile compounds such as water, residual solvents or added oils are evolved at relatively low temperatures. The elimination of such components depends on the gas pressure. At low pressures (vacuum), the corresponding mass loss step is shifted to lower temperatures, that is, vaporization is accelerated. The analysis of pyrolysis reactions in an inert atmosphere allows the content (from the step height) and possibly even the type of material to be determined.

The carbon black or carbon fiber content of a sample can be determined from the height of the combustion step after switching to an oxidative atmosphere. The residual filler, glass fiber or ash is determined from the residue. Small changes in the measurement curve due to buoyancy effects and gas flow rate can

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be corrected by subtracting a blank curve.

TGA measurements are often displayed as the first derivative of the TGA curve, the so-called DTG curve. Steps due to loss of mass in the TGA curve then appear as peaks in the DTG curve. The DTG curve corresponds to the rate of change of sample mass.

The temperature range of the decomposition steps is influenced to a certain extent by the ease with which the gaseous products are able to diffuse out of the sample. When reactive atmospheres are used, the efficiency of gas exchange at the surface of the sample is crucial. The effects of diffusion on the measurement can be reduced by using suitable crucibles (e. g. crucibles with low wall-heights such as the 30- μ L alumina crucible) and by suitable sample geometry (several small pieces or powder).

In TGA, the change in mass of the sample is measured very accurately. Unfortunately, however, the technique does not provide any information about the nature of the gaseous decomposition products evolved. The products can however be analyzed by coupling the TGA to a suitable gas analyzer (evolved gas analysis, EGA).

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Thermomechanical analysis measures the dimensional changes of a sample as it is heated. In this technique, the position or displacement of a probe resting on the surface of the sample with a certain force is continuously measured as a function of temperature or time. Figure 1.7 shows a typical TMA curve. The pressure exerted by the probe and the hardness of the sample determine whether the TMA experiment is in fact an expansion or a penetration measurement.

In the thermal expansion measurement, the probe exerts only a low pressure on the surface of the sample. The sample is heated linearly over the temperature range of interest. The linear coefficient of thermal expansion (CTE) is calculated directly from the measurement curve.

In a penetration experiment, the probe exerts a much greater pressure. The softening temperature can be directly measured when the sample is heated. Materials soften at the glass transition temperature or on melting.

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Fig. 1.7 Schematic TMA curve of a polymer under low compressive stress; 1 expansion below the glass transition; 2 glass transition (change of slope); 3 expansion above the glass transition; 4 plastic deformation

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If a periodically changing force is applied to the sample, the sample dimensions also change periodically. This measurement mode is called dynamic load TMA, DLTMA, and provides information on the viscoelastic behavior of polymers. The elastic modulus (Young's modulus) of the sample can be estimated from the amplitude and the sample thickness.

It is also possible to measure the bending behavior of hard samples with this technique using a bending accessory.

Special sample holders are available that allow the dimensional changes of fibers and films to be measured. A particular application is the measurement of the swelling behavior of materials in solvents. To do this, a sample is placed in a small container and the solvent of interest is added. The change in thickness of the sample due to solvent absorption is then continuously measured isothermally as a function of time.

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In dynamic mechanical analysis, a mechanical modulus is determined as a function of temperature, frequency and amplitude.

A periodically changing force (usually sinusoidal) applied to the sample creates a periodic stress in the sample. The sample reacts to this stress and the instrument measures the corresponding deformation behavior. The mechanical modulus, M , is determined from the stress and deformation. Depending on the

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type of stress applied, either the shear modulus, G (with shear stress) or the Young's modulus, E (with stretching or bending) is measured.

The sample does not always immediately react to the periodically changing stress—a certain time delay occurs that depends on the viscoelastic properties of the sample. This is the cause of the phase shift between the applied stress and the deformation. To take this phase shift into account, the dynamically measured modulus is described by a real part M' and an imaginary part M'' . The real part (storage modulus) describes the response of the sample in phase with the periodic stress. It is a measure of the (reversible) elasticity of the sample. The imaginary part (loss modulus) describes the component of the response that is phase-shifted by 90° . This is a measure of mechanical energy converted to heat (and therefore irreversibly lost). The tangent of the phase shift, $\tan \delta$, is also known as the loss factor and is a measure of the damping behavior of the material. The modulus and $\tan \delta$ depend on the temperature and the measuring frequency. At room temperature, rubbery materials show typical storage modulus values between 0.1 MPa and 10 MPa.

The curves in Figure 1.8 show the typical behavior of the storage and loss modulus of an amorphous and a semicrystalline polymer as a function of temperature.

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Fig. 1.8 Typical curves of the storage component (continuous line) and the loss component (dashed line) of the modulus as a function of the temperature for an amorphous polymer (blue) and a shock-cooled partially crystalline polymer (red): 1 secondary relaxation; 2 main relaxation; 3 rubbery plateau; 4 viscous flow; 5 cold crystallization; 6 melting

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At low temperatures the material is in a glassy state. The modulus is relatively high (about 2 GPa). In this state, a