Application Handbook Thermal Analysis

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This application handbook presents selected application examples. The experiments were conducted with the utmost care using the instruments specified in the description of each application at METTLER TOLEDO Thermal Analysis Lab in Switzerland. The results have been evaluated according to the current state of our knowledge.

This does not however absolve you from personally testing the suitability of the examples for your own methods, instruments and purposes. Since the transfer and use of an application is beyond our control, we cannot of course accept any responsibility.

When chemicals, solvents and gases are used, general safety rules and the instructions given by the manufacturer or supplier must be observed.

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Cyril Darribère

Preface

Over the past decades, the development and characterization of materials has become increasingly specialized due to the ever-increasing demands placed on quality. In addition, the widespread use of products such as blends, composites, new stabilizers and additives has resulted in relatively complex chemical systems. The great interest in material sciences has created a need for specific analytical tools to characterize the constituent inorganic and organic molecules. One particularly important development has been to combine well-established separation and identification techniques into one instrument system. This drastically increases the amount of specific information that can be obtained and shortens analysis times.

This booklet provides an insight into two such so-called hyphenated techniques, TGA-FTIR and TGA-MS. The first part of the booklet focuses on basic principles and describes the techniques. It also includes a practical section and an introduction to the interpretation of spectra. The second part discusses some 15 different applications performed in our applications laboratory using TGA-MS and/or TGA-FTIR. Two additional applications illustrate the use of the rather unusual combination of TMA with MS.

I would like to thank the Materials Characterization Section of the Market Support Group in Schwerzenbach, Switzerland for their helpful discussions, Mrs. Helga Judex for preparing the layout and Dudley May for checking the text.

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缩写和简称 Abbreviations and Acronyms	1
TGA-EGA 概述 Introduction to TGA-EGA	2
TGA-MS	4
TGA-FTIR	12
应用介绍 Introduction to the Applications	23
应用一览 Applications List	24
1 乙酰水杨酸分解 Decomposition of Acetylsalicylic Acid	25
2 BHET 的热裂解 Thermal Degradation of BHET	28
3 加热速率对 MS 响应的影响 Influence of the Heating Rate on MS Response	31
4 工业月桂醇的分解 Decomposition of Technical Lauryl Alcohol	34
5 药物中残留溶剂的检定 Detection of Residual Solvents in a Pharmaceutical Substance	37
6 一水草酸钙(演示品)分解 Decomposition of Calcium Oxalate Monohydrate (Tutorial)	39
7 样品量对 MS 灵敏度的影响 Influence of Sample Weight on MS Sensitivity	42
8 五水硫酸铜(演示品)分解 Decomposition of Copper Sulfate Pentahydrate (Tutorial)	44
9 填充有机物沸石的解吸附 Desorption of a Zeolite Filled with Organic Matter	47
10 PVC 粉末的分解 Pyrolysis of PVC Powder	50
11 氟化电缆线的研究 Investigation of Fluorinated Cable Wires	53
12 橡胶样品中水杨酸甲酯的检定 Detection of Methyl Salicylate in a Sample of Rubber	56
13 硅聚合物的裂解 Degradation of a Silicone Polymer	58
14 一种胺基树脂的固化和分解 Curing and Decomposition of an Amino Resin	61
15 借助热转变和分解鉴别 BR 和 NR 橡胶 Identification of BR and NR Rubbers by Using Thermal Transitions and Decomposition	65
16 印刷电路板的分层 Delamination of Printed Circuit Boards	69
17 用于包覆发泡剂的聚合物的膨胀和分解 Expansion and Decomposition of the Polymer Used to Encapsulate a Blowing Agent	71
参考文献 Suggested Literature	73

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Evolved gas analysis, EGA, is the term used to describe the analysis of gases or volatile components that are evolved from a sample undergoing thermal analysis.

Thermogravimetric analysis (TGA) is a quantitative technique used to measure the change in weight of a sample as a function of temperature (dynamically) or time (isothermally) under a controlled atmosphere. The simultaneously recorded SDTA signal (single differential thermal analysis) provides additional information on thermal effects such as evaporation, endothermic and exothermic chemical reactions, and phase transitions.

TGA alone, however, is not an identification technique. It cannot identify or characterize the nature of the gaseous products evolved during a thermogravimetric measurement. The combination of TGA with a mass spectrometer (MS) or a Fourier transform infrared (FTIR) spectrometer however makes this possible. The two techniques provide complimentary information and the data obtained can be correlated directly with the measured weight losses.

According to ICTAC (International Confederation of Thermal Analysis and Calorimetry), Evolved Gas Analysis (EGA) is defined as “a technique in which the nature and/or the amount of gas or vapor evolved from a sample is monitored against time or temperature while the temperature of the sample, in a specified atmosphere, is programmed (...). This technique involves a furnace and a gas analyzer or detector. Most commonly, the evolved gases are analyzed using mass spectrometry (MS) or Fourier transform infrared spectroscopy (FTIR)”.

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Time is the parameter common to the instruments making up the combination (i. e. TGA and MS, or TGA and FTIR). An

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interpretation of a measurement first of all requires a detailed study of the TGA curve. In addition to the quantification of the weight loss steps, the DTG curve, i. e. the first derivative of the TGA curve, allows the data obtained from the TGA to be correlated with that from the spectrometer. When the sample undergoes a change in weight, a peak maximum in the DTG curve (max. rate of change of weight) should correspond closely to a maximum in the signal intensity curve measured by the gas analyzer. Ideally, both MS and FTIR require the support of vapor phase spectral databases.

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


Mass spectrometry is an extremely sensitive method for the detection and identification of traces of gaseous substances. The technique can quantify atoms or molecules and provide chemical and structural information (functional groups and side chains) on the compounds analyzed.

In the mass spectrometers typically used for hyphenated techniques, sample molecules enter the mass spectrometer through an ion source where they are bombarded with a beam of high-energy electrons (usually about 70 eV). This energy is greater than the ionization potentials and bond strengths of organic substances. It is in fact sufficient to remove one (or more) electrons from molecules to form positively charged molecular ions. The excess energy also causes extensive fragmentation of the molecule with the formation of many different positively charged fragment ions according to complex pathways. The fragment ions formed are, however, directly related to the structure of the molecule under investigation, e. g.

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
The mass spectrometer separates the ions according to their mass-to-charge ratio m/z . With hyphenated techniques we are normally only interested in the singly charged molecule and fragment ions. For a number of reasons, quadrupole mass spectrometers have been widely used over the past twenty-five years in preference to conventional magnetic field sector analyzers in instrument combinations. A quadrupole mass analyzer consists of two pairs of parallel rod electrodes located between an ion source and a detector. The four rods are mounted opposite to each another and are connected in pairs (see diagram). Mass separation is achieved by the combination of a continuous electrostatic field and a high-frequency alternating field. The ions transmitted through the system are recorded at the SEM detector. Steadily increasing (scanning) the field strength causes ions of increasing m/z value to arrive at the detector

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and thereby produce a mass spectrum.

Mass spectrometers operate under high vacuum (typically 10^{-5} mbar). This is necessary to prevent ions colliding with residual gas molecules along their path from the ion source to the detector, which would prevent them from reaching the detector.

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
Schematic diagram of a quadrupole mass filter

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
Scan analog: The ion current is continuously measured as a function of the mass-to-charge ratio (i. e. for singly charged ions, the mass). The intensity and mass of the fragment ions provide information on the nature of the material under investigation.

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Scan paragraph: Only those masses are displayed whose signal intensity exceeds a particular value (i. e. that are above a certain threshold). The increment in m/z is therefore not continuous. The spectrum can be transformed into time-dependent intensities curves.

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
MID: Multiple Ion Detection (also known as Selected Ion Monitoring, SIM) enables you to continuously monitor the intensity of ions at defined m/z values, e. g. as shown in the following diagram at m/z 78 and 148.

In this technique the whole mass spectrum is not scanned, only a selected number of ions with particular m/z values are analyzed and the instrument rapidly switches between the selected ions until the measurements are completed. The amount of time (the so-called dwell time) that the analyzer remains at a given m/z value is thereby greatly increased. This in turn increases the fraction of these ions that reach the detector and results in a major improvement in sensitivity. Thus while subnanogram (10^{-10} to 10^{-9} g) quantities of compounds can normally be detected using full-scan quadrupole mass spectra, MID or SIM can in favorable cases lower the limit of detection to the picogram (10^{-12} to 10^{-11} g) range. MID is often applied to mixtures of compounds where sensitive quantitative or qualitative analysis is required.

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
A mass spectrum consists of a large number of signals of different intensity corresponding to the abundance of the molecular ion and the different fragment ions. It is therefore characteristic of the particular substance under investigation. Some important elements have more than one isotope with significant abundance (e. g. C, Si, Cl, Br). These isotope peaks are also recorded and often yield valuable additional information. Unfortunately, in general, the resolving power of the spectrometers used for hyphenated techniques is not sufficient to distinguish between ions of the same nominal m/z ratio (e. g. N^{2+} and CO^+ , both appear at m/z 28).

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The TGA/SDTA851[®] is coupled to the MS via a fused silica capillary enclosed in a heated transfer line. Since the MS operates at high vacuum, only about 1% of the effluent gas from the TGA is allowed to pass to the mass spectrometer (otherwise the vacuum would collapse). This small amount is however perfectly adequate because the sensitivity of the MS is extremely high. The remaining gas from the TGA is vented (to a collection tube). The transfer capillary is usually maintained at about 200°C to prevent condensation of the gaseous products. The pressure drop between the TGA (usually atmospheric pressure) and the MS (under high vacuum, $p < 10^{-5}$ mbar) causes the effluent from the TGA to be sucked into the MS.

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As mentioned above, a mass spectrum is often displayed as a bargraph, where the abscissa is the m/z ratio and the ordinate is the intensity of each ion. In addition, the signal intensities of selected fragment ions with particular m/z ratios (MID/SIM) can be also displayed as function of time (this aids correlation with a TGA or DTG curve). It is of course of the utmost importance that the masses assigned in the spectrum are correct—an error of just one mass unit would seriously confuse the interpretation.

The following guidelines are useful for the interpretation of mass spectra:

1. The abundance (intensity) of peaks often indicates whether the compound is aromatic (stable molecular ion, M^+) or aliphatic.
2. The fragment ion with the highest m/z generally corresponds to the ionized molecule, i. e. the molecular ion so that m/z is equal to M , the molecular weight of the compound. The ion peaks resulting from the loss of neutral fragments from the molecular ion peak are summarized in the following table:

从分子离子中失去中性碎片 Loss of neutral fragment from molecular ion	常见的中性碎片 Common neutral fragments
M-15	• CH ₃
M-16	• NH ₂
M-17	• OH
M-26	• CN
M-27	• HC=CH ₂
M-28	CO, C ₂ H ₄
M-29	• CH ₂ CH ₃ , • CHO
M-30	NO, NH ₂ CH ₂ •
M-31	• OCH ₃ , • CH ₂ OH
M-35	• Cl
M-41	CH ₂ =CHCH ₂ •
M-43	CH ₃ C • =O
M-45	• OCH ₂ CH ₃
M-49	• CH ₂ Cl
M-91	C ₆ H ₅ -C • H ₂

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3. If a compound contains only C, H, O, S and Si, the nominal m/z ratio of its molecular ion will be an even number; if a compound contains an even number of nitrogen atoms, the m/z of its molecular ion is an even number (nitrogen rule).

4. The relative intensities of isotopes (two atoms having the same atomic number but different atomic masses) can be very informative. For example, the relative abundance of ^{12}C compared to its isotope ^{13}C is 100:1.1. If a compound contains n carbon atoms and the molecular ion peak is normalized to 100, then the $M+1$ isotope peak is n times 1.1. It follows from this, that if no interferences occur, the maximum number of C atoms that an unknown molecule can contain can be calculated.

5. Isotopic effects are very helpful for identifying compounds containing chlorine and bromine. Chlorine has isotopes of ^{35}Cl and ^{37}Cl . Their relative natural abundances are 75% and 25% (3:1). This means that any compound containing one chlorine atom, e. g. chloromethane, shows a pair of peaks separated by two mass units with a 3:1 peak height ratio.

6. The m/z ratio of a peak in a spectrum does not exclusively identify a molecule or fragment ion, e. g. m/z 18. Water has a molecular weight of 18, but ionization can also produce ions that have the same nominal m/z ratio, e. g. the ammonium ion $[\text{NH}_4]^+$. Further information is required, e. g. in the latter case there should be peaks at m/z 14 and 15 corresponding to N^+ and $[\text{NH}]^+$ (atomic mass of nitrogen is 14).

7. The following pairs should be checked to avoid possible confusion in interpretation. This also means that careful consideration should be given to the choice of the purge gas:

m/z 17: OH and NH_3

m/z 28: N_2 and CO

m/z 29: $^{14}\text{N}^{15}\text{N}$ and C_2H_5 (alcohol)

m/z 36: ^{36}Ar and HCl

m/z 44: CO_2 and N_2O

8. The mass spectrum of an unknown compound consists of a number of peaks of different m/z ratios arising from the fragmentation of the initial molecule. The fragmentation processes can be classified as direct cleavage ($\text{ABCD}^+ \rightarrow \text{AB}^+ + \text{CD}$) or rearrangement ($\text{ABCD}^+ \rightarrow \text{AD}^+ + \text{B} = \text{C}$). When an