Application Handbook Thermal Analysis

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This application handbook presents selected application examples. The experiments were conducted with the utmost care using the instruments specified in the description of each application at METTLER TOLEDO Thermal Analysis Lab in Switzerland. The results have been evaluated according to the current state of our knowledge.

This does not however absolve you from personally testing the suitability of the examples for your own methods, instruments and purposes. Since the transfer and use of an application is beyond our control, we cannot of course accept any responsibility.

When chemicals, solvents and gases are used, general safety rules and the instructions given by the manufacturer or supplier must be observed.

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Jürgen E. K. Schawe 博士

Preface

This applications booklet provides an introduction to the thermal analysis of elastomers with a large number of practical examples. The main techniques used for sample measurement are thermogravimetric analysis (TGA), differential scanning calorimetry (DSC) and dynamic mechanical analysis (DMA). In particular cases, thermomechanical analysis (TMA) and combined techniques (TGA gas analysis) have also been employed.

In practice, thermal analysis is widely used for qualitative analyses. In fact however, the interpretation of the results and the final conclusions that can be drawn are greatly improved when quantitative measurements and evaluations are performed. We have therefore focused our attention on quantitative thermal analysis and outlined the basic principles for performing reliable measurements. The influence of experimental conditions on the measurement results is also covered.

Besides elementary considerations, a number of complex analyses and relationships are discussed.

I hope that the applications described in this book will find wide interest and stimulate new ideas both for experts and for newcomers to this rather complex but immensely interesting field.

I would like to thank Dr. Claus Wrana (Bayer AG, Leverkusen) in particular for making the samples available and for many discussions. I also thank Willi Weber (Huber and Suhner AG, Pfäffikon) for supplying the measurement data for Sections 4.3.1 and 4.7.1 and Christian Krebs (Uetendorf) for assistance in procuring samples.

Furthermore I would like to thank all my colleagues at METTLER TOLEDO, Schwerzenbach, Switzerland for their help and motivation; above all Dr. Markus Schubnell (for measurements on filled SBR and MS measurements), Cyril Darribère (for FTIR measurements), Dr. Rudolf Riesen (for TMA measurements and helpful discussions), Georg Widmann and Dr. Dudley May for proofreading the original German manuscript. I also owe my thanks to Dr. Dudley May for translating the German text into English and for many helpful suggestions.

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标题 Title	主题 Topics							方法 Methods				页码 Page
	玻璃化转变 Glass transition	结晶度/结晶 Crystallinity/crystallization	熔融 Melting	硫化 V vulcanization	组分 Composition	计算方法/实验 Evaluation/experimental Calculation/experiment	其他 Others	DSC/ADSC/IsoStep	TGA/TGA-EGA	TMA/DLTMA	DMA	
EPDM/SBR 共混物的 TGA 和 DSC 联合分析 Combined TGA and DSC analysis of an EPDM/SBR blend	•				•	•	•	•	•			213
未填充 SBR/NR 弹性体与温度关系的 DMA 测量 Temperature-dependent DMA measurements of an unfilled SBR/NR elastomer	•										•	218
未填充弹性体共混物的力学谱 Mechanical spectra of an unfilled elastomer blend	•						•				•	219
填充弹性体共混物与温度关系的 DMA 测量 Temperature-dependent DMA measurements of filled elastomer blends	•						•				•	221
填充弹性体共混物与频率关系的测量和主曲线绘制 Frequency-dependent measurements of a filled elastomer blend and master curve construction	•						•				•	223
升温速率对油含量测定的影响 The influence of the heating rate on oil content determination					•	•			•			226
含油和不含油 SBR 的降压(真空)TGA TGA of SBR with and without oil at reduced pressure (vacuum)					•	•			•			228
压力对 NR/SBR 共混物 TGA 的影响 Influence of pressure on the TGA of an NR/SBR blend					•		•		•			230
SBR 增塑的 DSC 测量 DSC measurements of the plasticization of SBR	•							•				232
三水合铝和氢氧化镁的 TGA TGA of aluminum trihydrate and magnesium hydroxide							•		•			233
EVA 中阻燃剂的 TGA TGA of flame retardants in EVA					•		•		•			235
阻燃物质的 DSC 测量 DSC measurements of flame-retarding substances							•	•				236
矿物油的 DSC 测量 DSC measurements of a mineral oil	•	•	•			•	•	•				237
软化测定的针入 TMA 测量 Determination of softening by penetration measurements by TMA	•						•			•		239
恢复行为的测量 Measurement of recovery behavior							•			•		241
溶胀行为的测量 Measurement of swelling behavior						•	•			•		242
密封圈适用温度的测定 Determination of the usable temperature range of sealing rings	•		•				•	•				244
氧化稳定性的等温 DSC 测量 Isothermal DSC measurements of oxidation stability						•	•	•				246

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应用一览表 Applications List	IV
1.热分析概论 Introduction To Thermal Analysis	1
1.1 差示扫描量热法 Differential Scanning Calorimetry (DSC)	1
1.1.1 常规 DSC Conventional DSC	1
1.1.2 温度调制 DSC Temperature-modulated DSC	2
1.2 热重分析 Thermo Gravimetric Analysis (TGA)	4
1.3 热机械分析 Thermo Mechanical Analysis (TMA)	6
1.4 动态热机械分析 Dynamic Mechanical Analysis (DMA)	7
1.5 与 TGA 的同步测量 SimuLTaneous Measurements with Tga	9
1.5.1 同步 DSC 和差热分析 Simultaneous DSC and differential thermal analysis (DTA, SDTA)	9
1.5.2 逸出气体分析 Evolved gas analysis (EGA)	11
1.5.2.1 TGA-MS	11
1.5.2.2 TGA-FTIR	12
2.弹性体的结构、性能和应用 Structure, Properties and Applications of Elastomers	14
2.1 聚合物的结构 Structure of Polymers	15
2.1.1 大分子中单体的排列 Arrangement of monomers in a macromolecule	15
2.1.2 大分子的结构 Structure of macromolecules	16
2.2 聚合物的物理结构 Physical Structure of Polymers	18
2.3 聚合物材料的分类 Classification of Polymeric Materials	20
2.4 弹性体产品 Elastomeric Products	22
2.4.1 聚合物的交联 Cross-linking of polymers	22
2.4.2 弹性体的生产 Manufacture of elastomers	23
2.4.3 弹性体的组成 Composition of elastomers	23
2.5 弹性体的分类 Classification of Elastomers	24
2.6 弹性体实例 Examples of Elastomers	25
2.7 弹性体分析 Analysis of Elastomers	28
2.7.1 热效应 Thermal effects	28
2.7.2 应用概述 Application overview	29
2.7.3 弹性体分析标准方法汇总 Summary of standard methods for elastomer analysis45	30
3.弹性体的基本热效应 Basic Thermal Effects of Elastomers	33
3.1 DSC 测试效应 Measurement Effects with DSC	33
3.1.1 玻璃化转变的计算方法 Evaluation of glass transitions	33
3.1.2 结晶和熔融 Crystallization and melting	42
3.1.3 硫化反应 Vulcanization reaction	50
3.1.4 借助 ADSC 用于曲线解析 The use of ADSC to aid curve interpretation93	60

3.2 效应的 TGA 测量 Measurement Effects with TGA	64
3.2.1 组分分析的基本原理 Basic principles of compositional analysis	64
3.2.2 实验条件的影响 Influence of experimental conditions	66
3.2.3 计算方法 Evaluation possibilities	68
3.3 TMA 测量 TMA Measurements	71
3.3.1 膨胀系数的测量 Measurement of the expansion coefficient	71
3.4 DMA 实验 DMA Experiments	73
3.4.1 复合力学性能 The complex mechanical properties	73
3.4.2 DMA 曲线的计算和表示 Evaluation and presentation of DMA curves	74
3.4.3 玻璃化转变和松弛 Glass transition and relaxation	82
4.热分析的应用 Applications Of Thermal Analysis	96
4.1 橡胶的分析 Analysis of Rubber	96
4.1.1 不同种类橡胶的 TGA 分析 TGA analysis of different types of rubber	96
4.1.2 玻璃化转变 Glass transition	106
4.1.3 结晶和熔融 Crystallization and melting	108
4.1.4 DSC 和 DMA 测量的比较 Comparison of DSC and DMA measurements	121
4.2 硫化 Vulcanization	123
4.2.1 硫化反应 Vulcanization reactions	123
4.2.2 SBR 硫化度对其热分析结果的影响 Influence of the degree of vulcanization on the TA results of SBR	133
4.2.3 硫化和未硫化硅橡胶的 DSC 比较 Comparison of vulcanized and unvulcanized silicone elastomers by DSC	147
4.3 碳黑的影响 The Influence of Carbon Black	153
4.3.1 碳黑的分析 Analysis of carbon black	154
4.3.2 测定碳黑对弹性体影响的 DSC 测量 DSC measurements to determine the effect of carbon black on elastomers	156
4.3.3 弹性体中碳黑的 TGA 分析 Analysis of carbon black in elastomers by TGA	158
4.3.4 含碳黑填料弹性体的 TMA 测量 TMA of elastomers with carbon black fillers	167
4.3.5 碳黑填充弹性体的 DMA 测量 Measurement of carbon black filled elastomers by DMA	169
4.4 以单一聚合物为基础的弹性体的实验 Experiments on Elastomers Based on a Single Type Of Polymer	171
4.4.1 TGA 分析 TGA Analysis	171
4.4.2 弹性体的气体分析鉴别 Identification of elastomers by means of gas analysis	180
4.4.3 DSC 测量 DSC Measurements	182
4.4.4 低分子量成分的熔融和结晶 Melting and crystallization of low molecular weight components	186
4.5 含若干个聚合物成分的弹性体 Elastomers with Several Polymer Components	191

4.5.1 TGA 测量 TGA Measurements	191
4.5.2 弹性体组分的 TGA/FTIR 鉴定 Identification of the polymer constituents by TGA-FTIR	203
4.5.3 DSC 测量 DSC measurements	205
4.5.4 复杂弹性体的分析实例 Example of a complex elastomer analysis	213
4.5.5 DMA 测量 DMA measurements	218
4.6 弹性体中油和增塑剂的分析 Analysis of Oils And Plasticizers In Elastomers	226
4.6.1 油含量的 TGA 测定 Determination of the oil content by TGA	226
4.6.2 增塑剂对玻璃化转变的影响 Influence of plasticizers on the glass transition	231
4.7 添加剂和助剂的分析 Analysis of Additives and Auxiliary Agents	233
4.7.1 阻燃剂 Flame retardants	233
4.7.2 增塑剂 Plasticizers	237
4.8 弹性体的进一步测量 Further Measurements of Elastomers	239
4.8.1 TMA 测量 TMA measurements	239
4.8.2 DSC 测量 DSC measurements	244
5.结果的总结 Summary of Results	248
参考文献 Literature	249

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Thermal analysis is the name given to a group of techniques that are used to measure the physical and chemical properties of materials as a function of temperature. In all these methods, the sample is subjected to a heating, cooling or isothermal temperature program.

The measurements can be performed in different atmospheres. Usually either an inert atmosphere (nitrogen, argon, helium) or an oxidative atmosphere (air, oxygen) is used. In some cases, the gases are switched from one atmosphere to another during the measurement. Another parameter sometimes selectively varied is the gas pressure.

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In DSC, the heat flow to and from the sample is measured relative to an inert reference material. The sample and reference materials are contained in separate but identical crucibles. This arrangement effectively compensates for thermal losses and the influence of the crucible.

DSC is used to measure thermal events such as physical transitions (glass transition, crystallization, melting, vaporization of volatile compounds) and chemical reactions. The information obtained characterizes the sample with regard to its thermal behavior and composition. In addition, with DSC it is also possible to measure quantities such as the heat capacity, glass transition temperature, melting temperature and degree of crystallinity.

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Conventional DSC employs a linear temperature program. The sample and reference material (or just an empty crucible) are either heated or cooled at a constant rate. Isothermal measurements can also be performed. Often several partial programs or segments are joined together to form a complete program. A schematic DSC curve is displayed in Figure 1.1. The change in the curve at the beginning of the measurement is due to the initial “startup deflection” (1). In this transient region, the conditions are suddenly changed from an isothermal to a linear heating mode. After the start up deflection, the sample is heated at the rate set for the measurement. The magnitude of the startup deflection depends on the heat capacity of the sample and the heating rate used. At a glass transition, the heat capacity of

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the sample increases and an endothermic step (2) is observed. A crystallization process (3) gives rise to an exothermic peak; the area of the peak corresponds to the enthalpy of crystallization. The melting of crystallites results in an endothermic peak (4). Chemical reactions give exothermic or endothermic peaks depending on the type of reaction involved. If volatile substances such as solvents are present in the sample, endothermic peaks (5) are observed due to vaporization; the sample consequently suffers a loss in mass. Further information on such peaks can be obtained by weighing the sample before and after the measurement and by using different types of crucibles. In contrast to open crucibles, the use of hermetically sealed crucibles prevents vaporization of the sample. Finally, at higher temperatures, decomposition begins (6). The type of purge gas used in the experiment often has a major influence on the reactions involved.

Transitions and reactions can be distinguished by measuring the same sample a second time after the sample has been cooled; chemical reactions are irreversible; crystalline materials however melt then crystallize again on cooling or on heating a second time.


图 1.1 图示 DSC 曲线:1 初始启动偏移;2 玻璃化转变;3 结晶;4 熔融;5 汽化;6 分解
Fig. 1.1 Schematic DSC curve: 1 initial startup deflection; 2 glass transition; 3 crystallization; 4 melting; 5 vaporization; 6 decomposition.

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Alternating DSC (ADSC) is a particular type of temperature modulated DSC. In contrast to conventional DSC, the temperature program is overlaid with a small periodic temperature change. The temperature program is characterized by the underlying heating rate,

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the temperature amplitude and the duration of the periodically changing temperature (Fig. 1. 2). With quasi-isothermal measurements, the underlying heating rate can also be zero.


图 1.2 典型 ADSC 温度程序: β_u 为基础升温速率; T_a 为温度振幅; t_p 为周期。
Fig. 1.2 Typical ADSC temperature program: β_u is the underlying heating rate, T_a the temperature amplitude, t_p period.


图 1.3 测得的 ADSC 曲线分离成基础和周期性信号成分。
Fig. 1.3 Separation of the measured ADSC curve into the underlying and the periodic signal components.

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As a result of temperature modulation, the measured heat flow changes periodically. This can be separated into two parts as shown in Figure 1. 3. Signal averaging yields the underlying signal (total heat flow), which corresponds to the conventional DSC curve. As additional information, one also obtains the periodic signal component. This is characterized by the amplitude of the heat flow and the phase shift between the actual (modulated) heating rate and the heat flow. The heat capacity (complex heat capacity) can then be determined from the amplitude and the phase shift. One advantage of this technique is