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Abstract

Preventive medical experiment has been an important practical course for students of clinical medicine at Tongji University. There is no officially published material and self-compiled teaching material has been used for a long time.

In order to standardize the experiment procedure and better use teaching material, we published this *Guide of Preventive Medical Experiment* on the basis of previous self-compiled book.

Guide of Preventive Medical Experiment aims to inspire the students' scientific thinking and creativity, help them apply theory to practice, and cultivate their ability to analyze and solve problems.

Use scope: mainly used by medical students in five-year, seven-year, eight-year programs and foreign students.

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Catalog

1. Determination of Formaldehyde in Indoor Air	5
2. Determination of Sulfur Dioxide Content in the Air	12
3. Disinfection of Drinking Water	18
4. Determination of Reductive-form Ascorbic Acid in Foods and Urine Load Test	25
5. Dietary Calculation and Evaluation	33
6. Discussion on Food Poisoning Cases	40
7. Neurobehavioral Function Test and Analysis	53
8. Diagnosis of Occupational Diseases	67
9. Discussion of Occupational Poisoning Case	77
10. Determination of Whole-Blood Acetylcholinesterase Activity	85

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Experiment 1

Determination of Formaldehyde in Indoor Air

[Objectives]

1. Master the significance of the determination of formaldehyde in indoor air.
2. Understand the principle and the determination method of formaldehyde in indoor air.

[Principle]

The formaldehyde in air reacts with phenol reagent to produce piperazine. In acidic solution, the piperazine is oxidized by high iron ion into blue-green compound. Quantitative analysis can be made according to the color.

[Materials]

Atmospheric sampling instrument;
10 ml Bubble absorption tube;
Spectrophotometer;
10 ml tube with a plug;
Pipets;
Aneroid barometer;
Absorption liquid;
Formaldehyde standard solution($1\text{ }\mu\text{g/ml}$);
1% Ferric sulfate amine solution.

[Procedures]

1. Add 5 ml of absorption liquid into a large bubble absorption tube, adjust gas velocity at 0.5 L/min, and sampling time is 20 min for 10L air.
2. Choose your own sampling sites. Record temperature and atmospheric pressure. After sampling, the sample should be analyzed within 24 h at room temperature.
3. Take 4 ml sample solution into 10 ml tube, and add 1 ml absorption solution, to make the total volume 5 ml. Mix well.
4. According to the following table, prepare the standard tubes series.

Table 1-1 Formaldehyde standard series

Tube number	0	1	2	3	4
Standard solution(ml)	0	0.2	0.4	0.8	1.6
Absorption solution(ml)	5.0	4.8	4.6	4.2	3.4
Formaldehyde content(μg)	0	0.2	0.4	0.8	1.6

5. Add 0.4 ml 1% ferric sulfate amine solution in the sample tube and standard tubes. Shake well and leave it for 15 min. With 1 cm thick cup, under the 630 nm, using blank solution to adjust zero, measure the optical density of each tube solution. Make standard curve first, and then the sample tube content is obtained from the standard curve.

6. Calculation

- (1) Convert the sample volume into sample volume in standard state.

$$V_0 = V_1 \times \frac{273 \times P}{(273 + t) \times 760}$$

Where,

V_1 ——sample volume(L);

P ——atmospheric pressure(mmHg);

T ——temperature($^{\circ}\text{C}$).

- (2) Calculate the formaldehyde concentration in the air

$$\text{Air formaldehyde concentration}(\text{mg}/\text{m}^3) = \frac{C}{V_0} \times \frac{5}{4}$$

Where,

C —sample content obtained from the standard curve (μg) ;

V_0 —sample volume in standard state(L).

[Notes]

1. Connect absorption tube with machine correctly.
2. Gas velocity should be stabilized at 0.5 L/min.
3. After adding reagent, mix well. The placing time is not less than 15 min.
4. Use the same color cup determinate the same set of samples and standard series.

[Discussions]

1. What is the main source of formaldehyde in indoor air?
2. What is the harm of formaldehyde on human health?

[Appendix]

The standard of formaldehyde in indoor air: 0.1 mg/m³ (average value of 1 hour).

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