APPLICATION HANDBOOK THERMAL ANALYSIS

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Jürgen de Buhr Georg Widmann / 著
..... 陆立明 / 译

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PHARMACEUTICALS & FOOD

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药物和食品 Pharmaceuticals & Foods

Jürgen de Buhr, Georg Widmann 著

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This application handbook presents selected application examples. The experiments were conducted with the utmost care using the instruments specified in the description of each application at METTLER TOLEDO Thermal Analysis Lab in Switzerland. The results have been evaluated according to the current state of our knowledge.

This does not however absolve you from personally testing the suitability of the examples for your own methods, instruments and purposes. Since the transfer and use of an application is beyond our control, we cannot of course accept any responsibility.

When chemicals, solvents and gases are used, general safety rules and the instructions given by the manufacturer or supplier must be observed.

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1 热分析概论 Introduction to Thermal Analysis	1
1.1 差示扫描量热法(DSC) Differential Scanning Calorimetry	1
1.1.1 常规 DSC Conventional DSC	1
1.1.2 温度调制 DSC(MTDSC) Temperature-modulated DSC	3
1.1.2.1 ADSC	3
1.1.2.2 IsoStep	4
1.1.2.3 TOPEM™	5
1.2 热重分析(TGA) Thermogravimetric Analysis	6
1.3 热机械分析(TMA) Thermomechanical analysis	8
1.4 动态热机械分析(DMA) Dynamic Mechanical Analysis	9
1.5 与 TGA 的同步测量 Simultaneous measurements with TGA	10
1.5.1 同步 DSC 和差热分析(DTA, SDTA) Simultaneous DSC and Differential thermal analysis ..	10
1.5.2 逸出气体分析(EGA) Evolved gas analysis	11
1.5.2.1 TGA/MS	11
1.5.2.2 TGA/FTIR	12
2 热分析在医药工业的应用 Applications of Thermal Analysis in the Pharmaceutical Industry	13
2.1 热分析药物应用一览表 Application Overview Pharmaceuticals	13
2.2 制药工业评说 Some Comments on the Pharmaceutical Industry	13
2.3 热分析在药物上的应用 Applications of Thermal Analysis in Pharmaceuticals	14
2.3.1 多晶型 Polymorphism	15
2.3.2 假多晶型 Pseudopolymorphism	15
2.3.3 相图 Phase diagrams	15
2.3.4 稳定性 Stability	16
2.3.5 相互作用 Interactions	16
2.3.6 纯度测定 Purity determination	17
2.3.7 包装材料 Packaging materials	17
2.3.8 工艺优化 Process optimization	17
2.3.9 校准和系统效应 Calibration, systematic effects	18
2.3.10 一些重要概念和缩写 Some important Concepts, Abbreviations and Acronyms	18
3 热分析的典型应用 Typical TA Applications of Pharmaceuticals	21
3.1 DSC 温度和热流量的校准 DSC Calibration, Temperature and Heat flow	21
3.2 与升温速率无关的 DSC 校准 DSC Calibration, Heating Rate Independence	22
3.3 升温速率对丁基羟基茴香醚多晶型检测的影响 Influence of Heating Rate on the Detection of Polymorphism, Butylated Hydroxyanisole	24
3.4 降温速率对蔗糖溶液结晶行为的影响 Influences of Cooling rate on Crystallization Behavior, Saccharose Solutions	25

3.5	升温速率对水包油乳膏水分含量测定的影响	Influences of Heating Rate on Moisture Content Determination, an O/W Cream	27
3.6	升温速率对美托拉胺分解的影响	Influences of Heating Rate on Decomposition, Metolazone	28
3.7	坩埚对一水葡萄糖失水的影响	Influence of the Pan on Dehydration, Glucose Monohydrate	29
3.8	丁基羟基茴香醚的样品制备	Sample Preparation, Butylated Hydroxyanisole	32
3.9	二丁基羟基甲苯试样量的影响	Influence of the Sample Weight, Butylated Hydroxytoluene	33
3.10	样品贮存和吸湿效应	Sample Storage and Hygroscopic Effects	34
3.11	油的氧化稳定性	Oxidation Stability of Oils	36
3.12	香草醛熔融行为的表征	Characterization of the Melting Behavior, Vanillin	37
3.13	胆固醇十四烷酸酯的相转变	Phase Changes, Cholesteryl Myristate	39
3.14	根据熔融行为对聚乙烯醇的鉴别	Identification Based on Melting Behavior, Polyethylene Glycol	40
3.15	糖溶液水的熔点降低	Melting Point Depression of Water, Sugar Solutions	42
3.16	油包水乳膏的 DSC“指纹”	DSC ‘Fingerprint’, O/W Cream	44
3.17	D,L 丙交酯-乙交酯共聚物的玻璃化转变	Glass Transition, Poly (D, L-lactide)-Co-Glycolide (DLPLGGLU)	45
3.18	羟丙基甲基纤维素邻苯二甲酸酯(HPMC-PH)的玻璃化转变和水分含量	Glass Transition and Moisture Content, Hydroxypropoxymethylcellulose Phthalate (HPMC-PH)	46
3.19	聚乙烯薄膜的质量控制	Quality Control, PE Films	49
3.20	氢化可的松的分解	Decomposition, Hydrocortisone	50
3.21	甲磺酸双氢麦角胺熔点处的分解	Decomposition at the Melting Point, Dihydroergotamine Mesylate	52
3.22	阿斯巴甜的熔融和分解	Melting Behavior and Decomposition, Aspartame	53
3.23	丙二酸的完全分解	Total Decomposition, Malonic Acid	55
3.24	乙酰水杨酸分解的动力学分析	Kinetic Analysis of Decomposition, Acetylsalicylic Acid	57
3.25	茶碱的水合稳定性	Hydrate Stability, Theophylline	59
3.26	淀粉/羟甲基纤维素钠(羧甲基淀粉钠)的水分	Moisture, Starch/NaCMC (Primojel)	61
3.27	三棕榈精的多晶型	Polymorphism, Tripalmitin	63
3.28	甲苯磺丁脲的多晶型	Polymorphism, Tolbutamide	64
3.29	退火处理丁基羟基茴香醚多晶型	Polymorphic Modifications by Annealing, Butylated Hydroxyanisole	66
3.30	硬脂酸镁的 DSC“指纹”	DSC ‘Fingerprint’, Magnesium Stearate	67
3.31	左旋聚丙交酯的多晶型	Polymorphism, L-Polylactide	69
3.32	磺胺吡啶的多晶型	Polymorphism, Sulfapyridine	70
3.33	一水葡萄糖的假多晶型	Pseudopolymorphism, Glucose Monohydrate	72
3.34	布洛芬(异丁苯丙酸)的光学纯度	Optical Purity, Ibuprofen	74

3.35	对羟基苯甲酸及其酯的纯度测定(DSC 法和 HPLC 法) Purity using DSC and HPLC, 4-Hydroxybenzoic Acid and its Esters	76
3.36	非那西汀+对氨基苯甲酸纯度测定 Purity Determination, Phenacetin + 4-Aminobenzoic Acid	78
3.37	胆固醇的纯度和重结晶 Purity and Recrystallization, Cholesterol	80
3.38	甲苯磺丁脲和聚乙二醇 6000 的相图 Phase Diagram, Tolbutamide and PEG 6000	82
3.39	对羟基苯甲酸甲酯和对羟基苯甲酸的共熔体组成 Eutectic Composition, Methyl-4- Hydroxybenzoate and 4-Hydroxybenzoic Acid	84
3.40	药物活性物质的 TGA-MS 溶剂检测 Solvent Detection by means of TGA-MS, Pharmaceutically Active Substance	86
3.41	不同水分含量的油包水乳膏的定量分析 Quantification, O/W Creams with Different Water Content	88
3.42	一水茶碱的定量分析 Quantification, Theophylline Monohydrate	90
3.43	Alcacyl 中活性物质的测定 Determination of an Active Substance, Alcacyl	92
4	热分析在食品工业的应用 Applications of Thermal Analysis in the Food Industry	95
4.1	热分析食品应用一览表 Application Overview Food	95
4.2	食品工业与热分析 Food Industry and Thermal Analysis	95
4.2.1	食品工艺中的反应和相 Reactions and Phases in Food Technology	95
4.2.2	食品中主要成分 DSC 检测一览表 List of DSC Investigations of the Main Components in Foods	96
4.2.3	蛋白质 Proteins	97
4.2.4	碳水化合物 Carbohydrates	99
4.2.5	脂肪和油 Fats and Oils	100
4.2.6	食品包装材料——塑料薄膜 Food Packagings-Plastic Films	102
5	热分析食品的典型应用 Typical TA Applications of Food	104
5.1	植物蛋白质的变性 Denaturation of Vegetable Proteins	104
5.2	鸡蛋蛋白质的变性 Egg Protein Denaturation	106
5.3	鸡蛋蛋清热处理的影响 Influence of Thermal Treatment of Egg White	108
5.4	鸡蛋贮存时间的影响 Influence of Egg Storage Time	109
5.5	pH 对牛血红蛋白的影响 Influence of pH on Bovine Hemoglobin	111
5.6	肉类的 DSC DSC of Meat	113
5.7	淀粉的凝胶化 Gelatinization of Starch	114
5.8	水中淀粉含量对溶胀的影响 Influence of the Starch Content on Swelling in Water	115
5.9	无定形糖的 ADSC(调制 DSC) ADSC of Amorphous Sugar	117
5.10	糖和淀粉的 TGA TGA of Sugar and Starch	119
5.11	意大利通心粉的动态负载 TMA Dynamic Load TMA of Pasta	121
5.12	巧克力的熔融 Melting of Chocolate	122
5.13	可可脂的热表征 Thermal Characterization of Cocoa Butter	124
5.14	熔融行为和氢化作用 Melting Behavior and Hydrogenation	128

5.15	植物油的结晶	Crystallization of Vegetable Oils	129
5.16	棕榈油的液相含量和滴点	Liquid Fraction and Dropping Point of Palm Oils	131
5.17	植物脂肪的氧化	Oxidation of Vegetable Fats	133
5.18	乙醇/水混合物	Ethanol/Water Mixtures	134
5.19	塑料薄膜的鉴别	Identification of Plastic Films	136

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Thermal analysis is the name given to a group of techniques used to measure the physical and chemical properties of materials as a function of temperature. In all these methods, the sample is subjected to a heating, cooling or isothermal temperature program.

The measurements can be performed in different atmospheres. Usually either an inert atmosphere (nitrogen, argon, helium) or an oxidative atmosphere (air, oxygen) is used. In some cases, the gases are switched from one atmosphere to another during the measurement. Another parameter sometimes selectively varied is the gas pressure.

DSC can also be used in combination with instruments that allow the sample to be simultaneously observed (DSC microscopy) or exposed to light of different wavelengths (photocalorimetry).

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In DSC, the heat flow to and from the sample is measured. DSC can be used to investigate thermal events such as physical transitions (the glass transition, crystallization, melting, and the vaporization of volatile compounds) and chemical reactions. The information obtained characterizes the sample with regard to its thermal behavior and composition. In addition, properties such as the heat capacity, glass transition temperature, melting temperature, heat and extent of reaction can also be determined.

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Conventional DSC employs a linear temperature program. The sample and reference material (or just an empty crucible) are heated or cooled at a linear rate, or in some cases, held at a constant temperature (i. e. isothermally). Often several partial programs or so-called segments are joined together to form a complete temperature program. A typical DSC curve is shown schematically in Figure 1.1. The change in the curve at the beginning of the measurement is due to the initial “startup deflection” (1). In this transient region, the conditions suddenly change from an isothermal mode to a linear heating mode. The magnitude of the startup deflection depends on the heat capacity of the sample and the heating rate. If volatile substances such as solvents are present in the sample, an endothermic peak (2) is

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observed due to the vaporization; the sample loses mass. Further information on such peaks can be obtained by weighing the sample before and after the measurement and by using different types of crucibles. In contrast to open crucibles, hermetically sealed crucibles prevent vaporization of the sample. At the part of DSC curve with no thermal effects (3), the heat capacity of the sample increases normally linearly and therefore a straight line called “baseline” is observed. The melting produces endothermic peak (4). Finally, at higher temperatures, decomposition begins (5). The type of purge gas used in the experiment often has an influence on the reactions that occur, especially at high temperatures.

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Figure 1.1 Schematic DSC curve

1 initial startup deflection; 2 evaporation of moisture; 3 part of DSC curve with no thermal effects, i. e., baseline; 4 melting peak; 5 beginning of decomposition.

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Physical transitions and chemical reactions can be differentiated by cooling the sample and measuring it again—chemical reactions are irreversible whereas crystalline materials melt then crystallize again on cooling or on heating a second time. Glass transitions are also reversible but not the enthalpy relaxation often observed in the first heating measurement of a glass transition.

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Alternating DSC (ADSC) is a particular type of temperature-modulated DSC (TMDSC). In contrast to conventional DSC, the linear temperature program is overlaid with a small periodic temperature change. The temperature program is characterized by the underlying heating rate, the temperature amplitude and the duration of the periodically changing temperature (Fig. 1.2). With quasi-isothermal measurements, the underlying heating rate can also be zero.

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
Figure 1.2 Typical ADSC temperature program

β_u is the underlying heating rate, A_T the temperature amplitude,

t_p period. The angular frequency ω is defined as $2\pi/P$ where P denotes the period of the sine wave.

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Figure 1.3 Separation of the measured ADSC curve into the underlying and the periodic signal components

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As a result of temperature modulation, the measured heat flow changes periodically. This can be separated into two parts as shown in Figure 1.3. Signal averaging yields the underlying signal (total heat flow), which corresponds to the conventional DSC curve. As additional information, one also obtains the periodic signal component. The reversing heat flow corresponds to the heat flow component that is able to follow the heating rate change directly and is computed from the in-phase heat capacity. The difference between the total heat flow and the reversing heat flow yields the non-reversing heat flow. One advantage of this technique is that it allows processes that occur simultaneously to be separated. For example, the change in heat capacity during a chemical reaction can be measured directly.

The evaluation of the ADSC curves is based on Fourier analysis. The modulus of the complex heat capacity c_p^* is calculated using the equation

$$|c_p^*| = \frac{A_\phi}{A_\beta} \cdot \frac{1}{m}$$

where A_ϕ and A_β denote the amplitudes of the modulated heat flow and heating rate, and m the sample mass. The phase angle between the ADSC heat flow signal and the heating rate is used to calculate the in-phase c_p .

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Fig. 1.4 IsoStep temperature program consisting of different isothermal and dynamic segments

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The isothermal segments allow the isothermal drift of the dynamic segments to be corrected. This results in better heat capacity accuracy. The isothermal step may also contain kinetic information, for example of a chemical reaction. Heat capacity determinations can be made using a sapphire reference sample, and kinetic effects can be separated from changes in heat capacity.

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TOPEM™ is an advanced temperature-modulated DSC technique that is based on the full mathematical analysis of the response of a DSC (both the apparatus and the sample) to a stochastically modulated underlying temperature program (Fig. 1.5). Due to the randomly distributed temperature pulses, the system is subjected to temperature oscillations over a wide frequency range and not just at one single frequency (e. g. , ADSC). An analysis of the correlation of the oscillating input signal (heating rate) and the response signal (heat flow) provides much more information than can be obtained using conventional temperature-modulated DSC. Not only can reversing and non-reversing effects be separated, but the quasi-static heat capacity of the sample is also measured and frequency-dependent heat capacity values are determined. This can be used to distinguish between frequency-dependent relaxation effects (e. g. glass transitions) and frequency-independent effects (e. g. chemical reactions) in one single measurement.

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Fig 1.5 Temperature curve of the furnace set value (black line) in a TOPEMTM method in which the furnace temperature (gray curve) generates a heating rate that fluctuates around a mean value.
The heat flow to the sample also fluctuates irregularly as shown in the diagram on the right.

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When a sample is heated, it often begins to lose mass. This loss of mass can result from vaporization or from a chemical reaction in which gaseous products are formed and evolved from the sample. If the purge gas atmosphere is not inert, the sample can also react with the gas. In some cases, the sample mass may also increase, e. g. in an oxidation reaction if the product formed is a solid.

In thermogravimetric analysis (TGA), the change in mass of a sample is measured as a function of temperature or time.

TGA provides information on the properties of the sample and its composition. If the sample decomposes as a result of a chemical reaction, the mass of the sample often changes in a stepwise fashion. The temperature at which the step occurs characterizes the stability of the sample material in the atmosphere used.

Figure 1.6 shows a typical TGA curve. The composition of a material can be determined by analyzing the temperatures and the heights of the individual mass steps.

Volatile compounds such as water, residual solvents or added oils are evolved at relatively low temperatures. The elimination of such components depends on the gas pressure. At low pressures (vacuum), the corresponding mass loss step is shifted to lower temperatures, that is, vaporization is accelerated. ? Loss of water of crystallization occurs usually above the temperature at which volatile compounds evaporate (2).

At higher temperature, the main component decompose and produces a rather high step (3). The gases used can be air or

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oxygen, or inert gases such as nitrogen. The analysis of pyrolysis reactions in an inert atmosphere allows the content (from the step height) and possibly even the type of material to be determined.

The filler or ash is determined from the residue. Small changes in the measurement curve due to buoyancy effects and gas flow rate can be corrected by subtracting a blank curve.

TGA measurements are often displayed as the first derivative of the TGA curve, the so-called DTG curve. Steps due to loss of mass in the TGA curve then appear as peaks in the DTG curve. The DTG curve corresponds to the rate of change of sample mass.

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Figure 1.6 Schematic TGA curve

1 loss of mass due to the vaporization of volatile components (moisture, solvents); 2 loss of water of crystallization; 3 decomposition; 4 residue (ash, fillers, carbon black or soot formed during decomposition in an inert atmosphere).

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The temperature range of the decomposition steps is influenced to a certain extent by the ease with which the gaseous products are able to diffuse out of the sample. When reactive atmospheres are used, the efficiency of gas exchange at the surface of the sample is crucial. The effects of diffusion on the measurement can be reduced by using suitable crucibles (e. g. crucibles with low wall-heights such as the 30- μL alumina crucible) and by suitable sample geometry (several small pieces or powder).

In TGA, the change in mass of the sample is measured very accurately. Unfortunately, however, the technique does not provide any information about the nature of the gaseous

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decomposition products evolved. The products can however be analyzed by coupling the TGA to a suitable gas analyzer (evolved gas analysis, EGA).

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Thermomechanical analysis measures the dimensional changes of a sample as it is heated. In this technique, the position or displacement of a probe resting on the surface of the sample with a certain force is continuously measured as a function of temperature or time. Figure 1.7 shows a typical TMA curve. The pressure exerted by the probe and the hardness of the sample determine whether the TMA experiment is in fact an expansion or a penetration measurement.

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In the thermal expansion measurement, the probe exerts only a low pressure on the surface of the sample. The sample is heated linearly over the temperature range of interest. The linear coefficient of thermal expansion (CTE) is calculated directly from the measurement curve.

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In a penetration experiment, the probe exerts a much greater pressure. The softening temperature can be directly measured when the sample is heated. Materials soften at the glass transition temperature or on melting.

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If a periodically changing force is applied to the sample, the sample dimensions also change periodically. This measurement mode is called dynamic load TMA, DLTMA. The elastic modulus (Young's modulus) of the sample can be estimated from the amplitude and the sample thickness.

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Figure 1.7 Schematic TMA curve

1 expansion below the glass transition; 2 softening (plastic deformation).