

**PROGRESS  
IN MEDICAL  
LABORATORY  
TECHNIQUE—3**

**F. J. BAKER**

PROGRESS IN  
MEDICAL LABORATORY  
TECHNIQUE

3

*Edited by*

F. J. BAKER

F.I.M.L.T., F.I.S.T., F.R.M.S.

*Chief Technician, Department of Pathology, Brompton  
Hospital; Associate Lecturer in Bacteriology,  
Paddington Technical College, London*

LONDON  
BUTTERWORTHS

1964

ENGLAND: BUTTERWORTH & CO. (PUBLISHERS) LTD.  
LONDON: 88 Kingsway, W.C.2

AUSTRALIA: BUTTERWORTH & CO. (AUSTRALIA) LTD.  
SYDNEY: 6-8 O'Connell Street  
MELBOURNE: 473 Bourke Street  
BRISBANE: 240 Queen Street

CANADA: BUTTERWORTH & CO. (CANADA) LTD.  
TORONTO: 1367 Danforth Avenue, 6

NEW ZEALAND: BUTTERWORTH & CO. (NEW ZEALAND) LTD.  
WELLINGTON: 49/51 Ballance Street  
AUCKLAND: 35 High Street

SOUTH AFRICA: BUTTERWORTH & CO. (SOUTH AFRICA) LTD.  
DURBAN: 33/35 Beach Grove

U.S.A.: BUTTERWORTH INC.  
WASHINGTON, D.C.: 7235 Wisconsin Avenue, 14

©

The several contributors named on pages vii and viii  
1964

*Printed in Great Britain by  
William Clowes and Sons, Limited, London and Beccles*

**URICASE** Leo  
for specific and accurate  
U/V spectrophotometric  
determination of uric acid

---

**URICASE-**  
specific uric acid-splitting  
enzyme

---

**URICASE-**  
the most accurate  
technique available

---

(standard deviation  
of 2 per cent.)\*

**URICASE-**  
uncomplicated, reliable,  
reproducible test

- \* **no protein precipitation**
- \* **no dangerous reagents**
- \* **simple enzymic reaction**
- \* **straightforward U/V readings**
- \* **a direct measure of blood uric acid**

Supplied in boxes containing  
1 ampoule of 75 units or  
5 ampoules of 15 units

\* Scand. J. clin. Lab. Invest., (1954), 6: 40

A 1489



Leo Laboratories Limited


London

Copenhagen · Dublin

# CENTRIFUGES


*Designed*  
for the  
**BIOCHEMIST**

**VERY LARGE  
CAPACITY**  
Up to 6 litres



'MISTRAL 6L'

**HIGH SPEED**  
Up to 18,000 r.p.m.  
combined with  
large capacity



'HIGH-SPEED 18'

**ULTRA SPEED**  
Up to 50,000 r.p.m.  
Precise temperature  
indication and control




'SUPER-SPEED 50'



**MEASURING & SCIENTIFIC EQUIPMENT LIMITED**  
Spenser Street, London, S.W.1, England. Telephone: VICTORIA 5426

# PAPANICOLAOU STAINS



**SHARP  
DEFINITION**

**In the  
cytological  
diagnosis  
of genital and  
extragenital  
cancer**

Ortho offer the following  
Papanicolaou staining  
materials:

**EA-50.** 100cc. and 500 cc.  
bottles

**EA-65.** 100cc. and 500 cc.  
bottles

**OG-6.** 100cc. and 500 cc.  
bottles


**Harris Haematoxylin  
ORTHO modification**  
100 cc. and 500 cc. bottles

**Vaginal cannulae and  
aspirator bulbs**  
for vaginal and  
endocervical smears

**Cervical spatulae  
(Ayre-type)**  
for cervical smears



LITERATURE ON REQUEST



**"GURR'S"**®

ESTABLISHED 1915

## **BIOLOGICAL STAINING METHODS**

by **George T. Gurr**



The *seventh edition* of this popular and practical manual. Most comprehensive publication of its kind Price 6/-

**STAINS  
AND  
REAGENTS  
FOR  
MICROSCOPY**

**I  
N  
C  
L  
U  
D  
I  
N  
G**

Orcein  
Silver Proteinate  
Coomassie Blue  
Fast Cyanin  
Dyes for Diagnostic Staining Procedures  
Dyes for Electrophoresis  
Buffer Tablets  
Immersion Oils  
Mounting Media  
Embedding Waxes  
Plastic Staining Vessels  
Embedding Baths and Ovens  
Microtomes  
Autoclaves  
Microscope Lamps  
Slide Storage Equipment

●  
THE CRITICS'  
CHOICE

U.K. AGENTS for  
*DANOCHEMO VITAMIN ASSAY MEDIA*  
*SERVÅ ION-EXCHANGE CELLULOSES*

●  
Write for Lists  
and Literature

**GEORGE T. GURR LTD.**  
136/144 NEW KINGS ROAD, LONDON, S.W.6

---

---

**E. & S. LIVINGSTONE, LTD.**

---

---

**THE CONTROL OF DISEASES**

By Hugh Paul, M.D., D.P.H.

*Second Edition.*

548 pages.

112 tables.

**75s.****ILLUSTRATED PHYSIOLOGY**

By Ann B. McNaught, M.B., Ch.B., Ph.D., M.I.Biol., and Robin Callander, F.F.Ph.

295 pages.

350 illustrations.

**30s.****ELECTRON MICROSCOPY**

By Gilbert Causey, M.B., F.R.C.S.

247 pages.

165 illustrations.

**48s.****CELL AND TISSUE CULTURE**

By John Paul, M.B., Ch.B., Ph.D., M.R.C.P.Ed., F.R.S.Ed.

*Second Edition.*

324 pages.

129 illustrations.

**32s. 6d.****HISTOPATHOLOGY**

By Robertson F. Ogilvie, M.D.(Edin.), D.Sc., F.R.C.P.(Edin.), F.R.S.E.

*Sixth Edition.*

528 pages.

400 illustrations in full colour.

**75s.****FLUORESCENT PROTEIN TRACING**

Edited by R. C. Nairn, M.D., Ph.D.

296 pages.

90 illustrations.

**42s.****ANATOMICAL TECHNIQUES**

By David H. Tompsett, B.Sc., Ph.D.

256 pages.

97 illustrations.

**35s.****THE INTERFERENCE MICROSCOPE IN BIOLOGICAL RESEARCH**

By Arthur J. Hale, M.B., Ch.B., D.Sc., Ph.D.

128 pages.

112 illustrations.

**20s.****CALLING THE LABORATORY**

Edited by W. A. R. Thomson, M.D.

124 pages.

**12s. 6d.****STAPHYLOCOCCUS PYOGENES AND ITS RELATION TO DISEASE**

By S. D. Elek, M.D., D.Sc., Ph.D., D.P.H.

775 pages.

81 illustrations.

**84s.****THE U.F.A.W. HANDBOOK ON THE CARE AND MANAGEMENT OF LABORATORY ANIMALS**

Edited by Alastair N. Worden, M.A., B.Sc., M.R.C.V.S., F.R.I.C., and

W. Lane-Petter, M.A., M.B., B.Chir.

*Second Edition.*

972 pages.

339 illustrations.

**70s.**

---

---

**TEVIOT PLACE, EDINBURGH**

---

---



not quite the  
thing!



... a rare chap like me requires the special made-to-measure care and attention that only a "Visijar" specimen container can give.

**"PERSPEX"**

(An I.C.I. main distributor)

We can manufacture all your "Perspex" requirements.

**"Visijar"**

TRADE MARK


"Perspex" Specimen Boxes and Apparatus

**THE "VISIJAR" LABORATORIES LTD.**

**149 London Road, Croydon, Surrey**

**Phone GRO 8228/0**

**Reagents for  
coagulation tests**



**Deprothrombinised  
Ox Plasma**

(as used for OWRENS  
P&P method)

Price: £8 per 100 vials  
each containing  
2 ml Ox Plasma  
£4 per 50 vials each  
containing 2 ml  
Ox Plasma  
£2 per 25 vials each  
containing 2 ml  
Ox Plasma

**Rabbit Brain  
Thromboplastin  
(Acetone Dried)**

Stability proven.  
Results Reproducible.  
Price: £1.15.0d. per box  
of six vials—each vial  
sufficient for 5ml  
of Thromboplastin.

**Rabbit Brain  
Thromboplastin  
(Saline Extract—  
Phenolized)**

Price: £1.15.0d per box  
of 6x5 ml bottles.


**Bell & Alton  
Platelet Substitute  
(For Thromboplastin  
Generation Tests)**

Price: 18/- per box  
of 6x5 ml bottles


  
Stayne Laboratories Ltd  
High Wycombe, Bucks

STAYNE  
STAYNE  
STAYNE  
STAYNE  
STAYNE  
STAYNE

Stayne Laboratories Ltd  
Holtbottom Road, Sandy  
High Wycombe, Bucks  
Tel: High Wycombe 5721



*Ortho Anti-Human Serum—  
the broad spectrum reagent  
with the confidence  
of anti-gamma and  
anti-non-gamma activity  
blended in every drop.*




## Anti-Human Serum



Ortho Pharmaceutical Limited  
Diagnostics Division  
Saunderton · Buckinghamshire

# LEEC

## ANIMAL CAGE STEAM STERILISING CABINET



★ **Another new addition to the  
LEEC range of Laboratory Apparatus**

This steam sterilising cabinet has been specially developed by LEEC and is giving excellent service in some of the leading laboratories. The internal measurements are: 42" high x 32" wide x 26" deep - made of stainless steel, with stainless-steel bars at the base. The exterior of the cabinet is in zinc-coated mild-steel sheet finished in epoxy enamel. The doors are insulated and there is a removable, debris collecting tray at the bottom of the cabinet.

*Write to the makers for full details:*


**LABORATORY & ELECTRICAL ENGINEERING CO.**  
71-73 Goldsmith Street, Nottingham Tel. Nottm. 48331

---

# **LEEC**

## **'PRECISION'**

### **LABORATORY INCUBATOR**



... designed by specialist engineers to meet the exacting demands of modern laboratory practice. The LEEC 'Precision' Incubator is now available with Circulating Fan, or with Cooling Coil and Fan. In addition, LEEC make a self-contained Refrigerator Unit which can be used with the LEEC Incubator with Cooling Coil. The Unit operates from ambient temperature to 0°C.

**TEMPERATURE  
VARIATION ...  
... less than 0.5°C in  
entire cabinet!**

*Write to the makers for full details:*

**LABORATORY & ELECTRICAL ENGINEERING CO.**  
71-73 Goldsmith Street, Nottingham Tel. Nottm. 48331

---

**TEKLAB (MEDICAL LABORATORIES) Ltd.**



**DISPOSABLE BLOOD COLLECTION TUBES**

*now available:*

- PROJECT I. *E.D.T.A.*
- PROJECT II. *Prothrombin*
- PROJECT III. *Urea*
- PROJECT IV. *Electrolyte*
- PROJECT V. *Blood Sugar*
- PROJECT VI. *Lithium Heparin 2ml.*
- PROJECT VIa. *Lithium Heparin 10ml.*
- PROJECT VII *E.S.R.*
- PROJECT VIII. *Serum Analyses*

*Full details and samples upon application*

**8 Gay Bowers, Hockley, Essex**


**Tel. 293**



## SUPPLIERS & STOCKISTS OF ALL YOUR LABORATORY REQUIREMENTS



including the



### UNIVERSAL JUNIOR II KS refrigerated centrifuge (illustrated)

Capacity from 4 × 800 ml. cups or blood bottles to 112 haematocrit tubes. Cup temperature can be maintained at 0°C even at top speeds. A wide range of heads and accessories available.

### also the Griffin SEROLOGICAL CENTRIFUGE for rapid blood grouping

... by the anti-globulin technique (Coombs, Mourant and Race, Lancet Feb 23, 1946, p. xxx). With this centrifuge, which is self-balancing, blood samples can be sedimented in less than one minute.

### and, for example, the

- New Griffin-Grundy Incubators
- Griffin 'Temperostat' and Water Bath
- Biochem Absorptiometer
- Unicam Colorimeters and Spectrophotometers
- E.I.L. Blood pH Meters
- Griffin van Slyke-neill Apparatus

Write today for the Griffin Comprehensive Catalogue of Modern Laboratory Apparatus

## GRIFFIN & GEORGE LTD.

Ealing Road, Alperton, Wembley, Middx. PERivale 3344


LONDON · BRISTOL · BIRMINGHAM · MANCHESTER · NEWCASTLE-UPON-TYNE · GLASGOW



Assures greater precision in the daily standardisation of your photometer and the periodic construction of your calibration curve. Each hermetically sealed ampoule contains an analytically measured, exact amount of cyanmethaemoglobin of assured purity and stability.



Facilitates the rapid preparation of diluent for accurate haemoglobin determinations by the cyanmethaemoglobin method. Simply dissolve a stable, standardised, sealed-in-foil pellet in 250 c.c. of distilled water



**NEW EASE  
AND ACCURACY  
IN HAEMOGLOBIN  
DETERMINATIONS**



ORTHO PHARMACEUTICAL LTD · DIAGNOSTICS DIVISION · SAUNDERTON · BUCKINGHAMSHIRE